

A human iPSC-based platform to screen therapeutics for ALS using assays covering disease-relevant readouts.



Benoit Samson-Couterie, Carola van Berkel, Leo Smit, Jessica Koepke, Shushant Jain
Ncardia Services BV, the Netherlands, www.ncardia.com


Background

Amyotrophic lateral sclerosis (ALS) is a rare and progressive neurodegenerative disorder that primarily affects motor neurons—the nerve cells responsible for controlling voluntary muscle movement. Clinically, patients experience increasing difficulty with speaking, swallowing, and mobility, and in advanced stages, respiratory failure due to loss of breathing muscle function. At the molecular level, most ALS cases are characterized by the mis-localization and pathological aggregation of the RNA-binding protein TDP-43 from the nucleus to the cytoplasm, where it forms insoluble inclusions. This disruption of normal TDP-43 function contributes to widespread defects in RNA processing. Among these, aberrant splicing and reduced expression of *STMN2*—a protein critical for axonal maintenance and neuronal repair—have emerged as key disease-associated phenotypes and biomarkers of TDP-43 dysfunction. Here, we demonstrate how human-relevant cellular assays can model these hallmark features of ALS pathology, including TDP-43 mis-localization, aggregation, and *STMN2* dysregulation. These platforms enable mechanistic studies and provide quantitative, translational tools to evaluate whether candidate therapies can rescue multiple aspects of disease pathogenesis, thereby supporting more predictive drug discovery efforts. Ncardia developed relevant assays using human iPSC derived motor neurons (MNs) with a TDP-43 point mutation to quantify ALS phenotypes:


- Mis-localization of TDP-43 to the cytoplasm
- Aggregation of TDP-43
- Reduction of *STMN2* protein levels
- Mis-splicing of *STMN2*
- Altered electrophysiological properties
- Neurofilament-L secretion

Methods


Ncardia employed hiPSC-derived MNs with a CRISPR engineered TDP-43 mutation* and assessed their ALS relevant phenotypes. All readouts were quantified after 2-3 weeks in culture.



Mislocalization of TDP-43 to the cytoplasm and TDP-43 aggregation



Quantifiable reduction of *STMN2* mRNA and mis-splicing of *STMN2*



Altered electrophysiological properties of TDP-43 MNs




Figure 3. Electrophysiological properties of wildtype and TDP-43 mutant MNs in coculture with hiPSC-derived astrocytes. (A) Quantification of electrophysiological activity in wildtype and TDP-43 mutant MNs, including active electrodes, mean firing rate (Hz), bursts, and network bursts. (B-C) Representative MEA raster plots of wildtype (B) and TDP-43 mutant (C) MNs co-cultured with hiPSC-derived astrocytes. Black lines = spikes, blue boxes = bursts, pink boxes = network bursts; each row represents one electrode (16 total). Data shown as mean ± SEM; *p<0.05, **p<0.005, ***p<0.0005, ****p<0.00005.

Small molecule rescue of TDP-43 aggregation and electrophysiological deficits




Figure 4. MEA assessment of network activity in wildtype and TDP-43 mutant motor neurons (MNs). (A) Representative MEA raster plots of MNs co-cultured with hiPSC-derived astrocytes. Each row represents one electrode (16 total). Black ticks = spikes, blue boxes = bursts, and network bursts. Shown are wildtype, TDP-43 mutant, and rescue/treated conditions. Quantification of active electrodes, mean firing rate (Hz), bursts, and network bursts. Data are presented as mean ± SD; *p<0.05, **p<0.005, ***p<0.0005, ****p<0.00005. (B) Quantification of TDP-43 aggregation shown as % change relative to control. Mutant MNs display increased aggregation, which is reduced by treatment. Data are mean ± SEM; **p<0.01, ***p<0.001, ****p<0.0001 vs control.

Successful transduction of TDP-43 mutant MNs & astrocyte coculture subsequent decrease of TDP-43 aggregation




Figure 5. Transduction of TDP-43 mutant MNs in coculture with hiPSC-derived astrocytes and significant reduction of TDP-43 aggregation. (A) Representative images of TDP-43 mutant MNs in coculture with hiPSC-derived astrocytes 1 and 2 weeks after transduction with GFP carrying gene therapy. (B) Quantification of TDP-43 aggregation in TDP-43 mutant MNs in coculture with hiPSC-derived astrocytes treated with gene therapy, normalized to untreated cultures. Error bars represent mean ± SEM; *p<0.05, **p<0.005, ***p<0.0005, ****p<0.00005.

Conclusions

- We demonstrate quantifiable TDP-43 mislocalization and *STMN2* protein reduction by high-throughput high-content imaging and miniaturized assays. Furthermore, this model exhibits other relevant hallmarks of ALS, such as quantifiable TDP-43 aggregation, *STMN2* mis-splicing, neurofilament-L secretion and electrophysiological deficits as measured by MEA.
- We present a human in vitro model amenable for evaluation of new therapeutic candidates. The cultures were successfully transduced and when treated with a gene therapy, Ncardia was able to significantly reduce TDP-43 aggregation as compared to non-treated cultures.

Comprehensive physiologically-relevant platform, with hiPSC-based miniaturized assays, for screening of ALS therapeutics

Ncardia Services B.V.

Contact us at
info@ncardia.com or
scan the QR code to
download the poster

