

Development of robust iPSC-based α -Synuclein, Tau and TDP-43 aggregation models for drug discovery



Jessica Koepke, Carola van Berkel, Leo Smit, Benoit Samson-Couterie and Shushant Jain


Ncardia Services BV Leiden, The Netherlands - No conflicts of interest

Background

Proteinopathies are disorders driven by protein misfolding and aggregation, resulting in impaired neuronal function and progressive neurodegeneration. Notable examples include Parkinson's disease (PD), Alzheimer's disease (AD), and Amyotrophic Lateral Sclerosis (ALS). These conditions are highly complex and remain challenging to model accurately in both in vitro and in vivo systems. Consequently, many therapies that show promise in animal models fail in clinical trials, highlighting a significant translational gap in drug discovery.

Developing physiologically relevant human disease models is therefore critical to more reliably evaluate therapeutic candidates and improve clinical success rates. Human induced pluripotent stem cells (hiPSCs) offer a powerful solution, as they can differentiate into diverse cell types, preserve patient-specific genetic backgrounds, and recapitulate human pathophysiology while responding predictively to candidate drugs. These features make hiPSC-based systems well suited for drug discovery applications.

In this study, Ncardia established three hiPSC-derived neuronal in vitro assays to model the aggregation of α -synuclein, Tau, and TDP-43—key pathological hallmarks of PD, AD, and ALS.




Postulated model for protein aggregation mechanisms in neurodegenerative diseases

- Misfolding of specific characteristic disease related proteins is suggested to be linked to disease progression, resulting in aggregation and fibril formation of these proteins
- Once disease protein aggregates and loses its function or additionally shows a toxic gain of function.


Bahareh Eftekhari, Bradley T. Hyman, Susanne Wegmann. Structural studies on the mechanism of protein aggregation in age related neurodegenerative diseases. *Mechanisms of Ageing and Development*, Volume 156, 2016, Pages 1-13, ISSN 0047-6374, <https://doi.org/10.1016/j.mad.2016.03.001>.

Aggregation and phosphorylation of α -synuclein




Treatment with α -syn pre-formed fibrils (PFF) results in p- α -syn aggregates

Electrophysiological dysfunction in a model of PD




Treatment of Ncyte CNS with α -syn monomers or PFF disturbed Ncyte CNS electrical activity

Quantification of pTAU, soluble and insoluble Tau



Overexpression of TAU combined with chronic exposure to PFFs induces an Alzheimer's disease phenotype

Mis-localization and aggregation of TDP-43 and downregulation STMN2



Motor neurons chronically treated with a stressor show mis-localization of TDP-43 to the cytoplasm, aggregation of TDP-43, reduction of STMN2 and appearance of the truncated variant

Conclusions

- Neuronal co-cultures were used to quantify disease relevant phenotypes for α -synuclein or TAU aggregation as well as the formation intermediate phosphorylated species after treatment with α -synuclein and Tau recombinant preformed fibrils (PFFs).
- Stressor-treated mutant and wild type iPSC-derived motor neurons (hiPSC-MN) showed disease-specific mis-localization of TDP-43 to the cytoplasm, aggregation of TDP-43, reduction of STMN2 and appearance of the truncated variant.

- We have established a suite of robust, clinically relevant in vitro assays (Z -factor >0.5). These assays are fully automated and performed in a scalable format to support drug discovery at any stage.

- We successfully modelled and evaluated disease-linked phenotypes relevant to AD, PD and ALS, among other neurodegenerative disorders, using complementary assays. Altogether, offering the opportunity to gain a holistic understanding of the efficacy of therapies targeting aggregation

Scan to download the poster

Questions? Contact us at: support@ncardia.com www.ncardia.com

