

# A human iPSC-based platform to screen therapeutics for ALS using specific and robust phenotypic assays



Carola van Berkel, Leo Smit, Benoit Samson-Couterie, Jessica Koepke, Shushant Jain

Ncardia Services BV, the Netherlands, support@ncardia.com, www.ncardia.com

## Background


Amotrophic lateral sclerosis (ALS) is a rare neurological disease that primarily affects the nerve cells (motor neurons) responsible for controlling voluntary muscle movement. ALS is characterized by stiff muscles, muscle twitching, and gradually worsening weakness due to muscles decreasing in size. This results in difficulty speaking, swallowing, and eventually breathing. Most of the cases present mis-localization and pathological aggregates of TDP-43 in the cytoplasm. Furthermore, mis-splicing of *STMN2* and subsequent reduction of *STMN2* protein are also core phenotypes associated with ALS.

Ncardia developed relevant assays using human iPSC derived motor neurons (MNs) with a CRISPR engineered TDP-43 point mutation to quantify ALS phenotypes:

- Mis-localization of TDP-43 to the cytoplasm
- Aggregation of TDP-43
- Reduction of *STMN2* protein levels
- Mis-splicing of *STMN2*
- Altered electrophysiological properties
- Neurofilament-L secretion

## Methods

Ncardia employed hiPSC-derived MNs with a CRISPR engineered TDP-43 mutation\* and assessed their ALS relevant phenotypes. All readouts were quantified after 2-3 weeks in culture.



\* iCell® Motor Neurons, ALS TDP43, 01279 and iCell® Motor Neurons, 01279 from FUJIFILM Cellular Dynamics, Inc

## Mislocalization of TDP-43 to the cytoplasm and pathological TDP-43 aggregate formation




Figure 1. Mislocalization of TDP-43 to the cytoplasm and pathological aggregation of TDP-43 in hiPSC-derived TDP-43 mutant MNs treated with a stressor

(A) HCl images (40x) of vehicle-treated compared to stressor-treated TDP-43 mutant MNs. Immunoreactivity to DAPI in blue, TDP-43 in red and MAP2 in grey. Zoom-in of relevant structures in bottom-right. (B) Quantification of the nuclear/cytoplasmic ratio of TDP-43 intensity, normalized to vehicle-treated MNs. (C) Quantification of TDP-43 aggregation of stressor-treated TDP-43 mutant MNs and hiPSC-derived astrocytes, normalized to vehicle-treated cultures. Error bars represent mean  $\pm$ SEM, \*p<0.05, \*\*p<0.005, \*\*\*p<0.0005, \*\*\*\*p<0.00005.

## Quantifiable reduction of *STMN2* protein and mis-splicing of *STMN2* gene




Figure 2. Reduction of *STMN2* protein and mis-splicing of *STMN2* gene

(A) HCl images (40x) of vehicle-treated compared to stressor-treated TDP-43 mutant MNs. Immunoreactivity to DAPI in blue, *STMN2* in green and MAP2 in grey. Zoom-in of relevant structures in bottom-right. (B) Quantification of *STMN2*  $\Sigma$  intensity per  $\mu\text{m}^2$  in neurites (defined by MAP2), normalized to vehicle-treated MNs. (C) Quantification of truncated *STMN2* mRNA of stressor-treated TDP-43 mutant MNs, normalized to stressor-treated cultures. Error bars represent mean  $\pm$ SEM, \*p<0.05, \*\*p<0.005, \*\*\*p<0.0005, \*\*\*\*p<0.00005.

## Altered electrophysiological properties of TDP-43 MNs in complex coculture with astrocytes




Figure 3. Electrophysiological properties of wildtype and TDP-43 mutant MNs in coculture with hiPSC-derived astrocytes

(A) Representative raster plots of multi-electrode array (MEA) recordings of wildtype MNs in coculture with hiPSC-derived astrocytes. Black lines indicate spikes, blue boxes indicate bursts and network bursts are indicated by pink boxes. Each row represents 1 electrode (total 16 electrodes). (B) Representative raster plots of MEA recordings of TDP-43 mutant MNs in coculture with hiPSC-derived astrocytes. Black lines indicate spikes, blue boxes indicate bursts and network bursts are indicated by pink boxes. Each row represents 1 electrode (total 16 electrodes). (C) Quantification of electrophysiological activity of wildtype and TDP-43 mutant MNs. Including number of active electrodes, mean firing rate (Hz), number of bursts and network bursts. Error bars represent mean  $\pm$ SEM, \*p<0.05, \*\*p<0.005, \*\*\*p<0.0005, \*\*\*\*p<0.00005.

## Neurofilament-L release in TDP-43 mutant MNs




Figure 4. Quantification of neurofilament-L release of TDP-43 mutant MNs

Quantification of Neurofilament-L (NF-L) in supernatants of TDP-43 mutant MNs treated with a stressor normalized to vehicle-treated TDP-43 MNs. Error bars represent mean  $\pm$ SEM, \*p<0.05, \*\*p<0.005, \*\*\*p<0.0005, \*\*\*\*p<0.00005.

## Successful transduction of TDP-43 mutant MNs & astrocyte coculture with gene therapy and subsequent reduction of TDP-43 aggregation




Figure 5. Transduction of TDP-43 mutant MNs in coculture with hiPSC-derived astrocytes and significant reduction of TDP-43 aggregation

(A) Representative images of TDP-43 mutant MNs in coculture with hiPSC-derived astrocytes 1 and 2 weeks after transduction with GFP carrying gene therapy. (B) Quantification of TDP-43 aggregation in TDP-43 mutant MNs in coculture with hiPSC-derived astrocytes treated with gene therapy, normalized to untreated cultures. Error bars represent mean  $\pm$ SEM, \*p<0.05, \*\*p<0.005, \*\*\*p<0.0005, \*\*\*\*p<0.00005.

## Conclusions

- In this poster we present our hiPSC-derived cellular model with ALS-relevant phenotypes all of which we have miniaturized to 384-well format suitable for high-throughput compound screening.
- We demonstrate quantifiable TDP-43 mislocalization and *STMN2* protein reduction by high-throughput high-content imaging. Furthermore, this model exhibits quantifiable TDP-43 aggregation, *STMN2* mis-splicing, neurofilament-L secretion and electrophysiological deficits as measured by MEA.
- Finally, we demonstrate that this model is amenable for evaluation of new therapeutic candidates. The cultures were successfully transduced and when treated with a gene therapy, Ncardia was able to significantly reduce TDP-43 aggregation as compared to non-treated cultures.

Comprehensive physiologically-relevant platform, with hiPSC-based miniaturized assays, for screening of new ALS therapeutics!

Scan this QR code to download the poster

