

Evaluation of human induced pluripotent stem cell (hiPSC)-derived tri-culture as in vitro model for neuroinflammation

Isabel Onofre, Chiara Formica, Leo Smit, Joyson Dcosta and Shushant Jain


Ncardia BV Leiden, The Netherlands

No conflicts of interest

Background

The interplay between the immune and nervous systems is a critical aspect of various neurodegenerative diseases (NDD), but the development of physiologically relevant models mimicking the intercellular interactions remains a challenge. To overcome this, human induced pluripotent stem cell technology can be leveraged to build multicellular models and bring human pathophysiology into early drug discovery to develop new therapeutics targeting neuroinflammation. The development and validation of models of relevant biological disease processes, such as microglia-neuron communication provides insight on cellular interactions that play a role in recognizing apoptotic neurons and modulating neuronal activity which are crucial events in disease progression. Targeting these pathways in human models with a combination of readouts allows interrogation and evaluation of the ability of therapeutics on rescuing primary, secondary and tertiary neuro-pathological signatures.

In this study, Ncardia used the tri-culture approach to developed an *in vitro*, iPSC-derived tauopathy assay by inducing phosphorylation (phospho-Tau) and aggregation of Tau, with recombinant mutant Tau PFFs (pre-formed fibrils) and Tau seeds extracted from AD patient brain tissue. This approach enabled a multi-parametric readout of neuronal and glial phenotypes including activation of microglia and astrocytes in the tri-culture.



Cross-talk between neurodegeneration and neuroinflammation

- Neuronal homeostasis relies on astrocytes and microglia activation status, that is disturbed in pathogenic conditions as neurodegeneration
- Diseased neurons secrete factors that activate microglia and astrocytes, initiating a cascade of inflammatory triggers that induce further neurodegeneration and neuroinflammation

Adapted from: Ullah, F., Gamage, R., Sen, M.K. et al. The Effects of Modified Curcumin Preparations on Glial Morphology in Aging and Neuroinflammation. *Neurochem Res* 47, 813-824 (2022). <https://doi.org/10.1007/s11064-021-03499-4>






Figure A shows the experimental set-up, timeline and endpoints of Ncardia's neuroinflammation model. Assays in purple can be combined and multiplexed with assays in green

1. Tri-culture model: neurons, astrocytes and microglia




Ncardia's tri-culture model showed functional activation of microglia upon stimulation with LPS

2. Induction of neuroinflammation as a model for AD using TAU PFFs



Ncardia's tri-culture model showed functional activation of microglia upon stimulation with Tau PFF

3. Induction of neurodegeneration and neuroinflammation as a model for AD with human brain extracts



Ncardia's tri-culture model showed neurodegeneration and neuroinflammation phenotype upon treatment with AD Tau seeds

Conclusions

- We have established a physiological relevant human *in vitro* tri-culture model, composed by neurons, astrocytes and microglia suitable to study the effect of drug candidates on neuroinflammation and neurodegeneration.
- Microglia is activated upon stimulation with LPS and Tau PFF, increasing the release of pro-inflammatory cytokines, a hallmark of neuroinflammation. Neurons treated with Tau PFFs expressed and accumulated phospho-Tau, exhibiting a neurodegenerative phenotype.

- Upon the addition of Tau seeds extracted from AD brain tissue, it was observed increased levels of expression of phospho-Tau (pTau) in neurons and released Tau in supernatants. Microglia phagocytosed neurons expressing phospho-Tau and released higher levels of Tau and IL-6. Together, these observations support a neurodegenerative phenotype, typical of tauopathies in which secreted or engulfed phospho-Tau activates microglia initiating and contributing to the neuroinflammatory cascade.
- These models are available on-demand to help drug developers select the best candidates earlier in the drug discovery process, increasing confidence and reducing the need for animal models.



Scan to download the poster

Questions? Contact us at: support@ncardia.com
www.ncardia.com