

Development of robust iPSC-based α -Synuclein, Tau and TDP-43 aggregation models for drug discovery

Jessica Koepke, Isabel Onofre, Carola van Berkel, Tatiana Wust, Leo Smit, Joyson Dcosta, Chiara Formica, Kimberley Riegman, Benoit Samson-Couterie and Shushant Jain
Ncardia Services BV Leiden, The Netherlands

No conflicts of interest

Background

Proteinopathies are diseases caused by protein misfolding and self-aggregation, which leads to altered neuronal function and neurodegeneration. Examples of proteinopathies include Parkinson's Disease (PD), Alzheimer's Disease (AD) and Amyotrophic lateral sclerosis (ALS). Proteinopathies are complex diseases difficult to model in vitro and in vivo. Many therapies that showed promise in animal studies have failed in human clinical trials, emphasizing a translation gap in the drug discovery process. Developing human physiologically relevant disease models is of high importance to identify and validate drugs therapeutic potential with higher confidence of clinical success. Human induced pluripotent stem cells (hiPSCs) have the potential to be differentiated into any cell type, retain patient-specific genetic backgrounds, mimic clinically-relevant human (patho-)physiology and respond appropriately to candidate therapeutics. These characteristics make hiPSCs an excellent tool for drug discovery.

In this study, Ncardia developed three in vitro assays based on hiPSC-derived neurons to model the aggregation of α -synuclein, Tau and TDP-43 – key hallmarks for diseases like Parkinson's, Alzheimer's or Amyotrophic Lateral Sclerosis.





Figure A shows the experimental set-up, timeline and endpoints of Ncardia's AD, PD and ALS pathology model. Assays in purple and in green can be combined and multiplexed




Postulated model for protein aggregation mechanisms in neurodegenerative diseases
• Misfolding of specific characteristic disease related proteins is suggested to be linked to disease progression, resulting in aggregation and fibril formation of these proteins

• One disease protein aggregates and loses its function or additionally shows a toxic gain of function.


Bahareh Eftekharaezadeh, Bradley T. Hyman, Susanne Wegmann. Structural studies on the mechanism of protein aggregation in age related neurodegenerative diseases. *Mechanisms of Ageing and Development*, Volume 156, 2016, Pages 1-13, ISSN 0047-6374, <https://doi.org/10.1016/j.mad.2016.03.001>.

1. Aggregation of α -synuclein




Treatment with α -syn pre-formed fibrils (PFF) results in p α -syn aggregates

2. Electrophysiological dysfunction in a model of PD




Treatment of Ncyte CNS with α -syn monomers or pre-formed fibrils (PFF) disturbed Ncyte CNS electrical activity

3. Aggregation of Tau and pTau - Quantifying aggregation



Overexpression of TAU combined with chronic exposure to PFFs induced an Alzheimer's disease phenotype

4. Mis-localization and aggregation of TDP-43 and downregulation STMN2



Motor neurons chronically treated with a stressor show mis-localization of TDP-43 to the cytoplasm, aggregation of TDP-43, reduction of STMN2 and appearance of the truncated variant

Conclusions

• Neuronal co-cultures were used to quantify disease relevant phenotypes for α -synuclein or TAU aggregation as well as the formation intermediate phosphorylated species after treatment with α -synuclein and Tau recombinant preformed fibrils (PFFs) using a high content imaging platform. Stressor-treated mutant and wild type iPSC-derived motor neurons (hiPSC-MN) showed disease-specific mis-localization of TDP-43 to the cytoplasm, aggregation of TDP-43, reduction of STMN2 and appearance of the truncated variant.

• We have established a suite of robust, clinically relevant in vitro assays (Z-factor >0.5) for the aggregation of α -synuclein, Tau and TDP-43 using human iPSC-derived neuronal subtypes. These assays are performed in a scalable wellplate format and are fully automated to support drug developers at any stage of their discovery process.

• We successfully modelled and evaluated disease-linked phenotypes relevant to AD, PD and ALS, among other neurodegenerative disorders, using complimentary assays. Altogether, offering the opportunity to gain a holistic understanding of the efficacy of therapeutics targeting aggregation.



Scan to download the poster

Questions? Contact us at:
support@ncardia.com
www.ncardia.com