

Evaluation of human induced pluripotent stem cell (hiPSC)-derived tri-culture as in vitro model for Alzheimer’s Disease

L. SMIT. VAN BERKEL, L. SMIT, J. COTOVIO, P. VAN LOENEN, M. ARGENZIANO , S. JAIN;

Ncardia Services BV Leiden, The Netherlands, support@ncardia.com, www.ncardia.com




Background

The development of physiologically relevant models for Alzheimer’s Disease (AD) remains a challenge with an unfilled gap in translatable human-based platforms. Ncardia developed human tri-culture models composed of neurons, astrocytes and microglia derived from iPSCs, resembling physiological conditions that allow modulation of neuroinflammation and neurodegeneration in vitro in the context of AD.

Firstly, we identified relevant triggers as TAU species capable of inducing cellular responses specific to AD pathology in cultures of microglia.


Secondly, we established a tauopathy model in the tri-culture system, by inducing the phosphorylation (pTAU), misfolding and aggregation of TAU, using different recombinant mutant TAU (pre-formed fibrils) PFFs. This approach enabled a multi-parametric readout of neuronal and glial phenotypes including activation of microglia and astrocytes in the tri-culture, and release of pro-inflammatory cytokines in the supernatants.

The development and validation of models of relevant biological disease processes, such as microglia-neuron communication, provides insight on cellular interactions. Modeling these cellular interactions play a role in recognizing apoptotic neurons and modulating neuronal activity, which are crucial events in disease progression. Targeting these pathways in human models with a combination of clinically-relevant readouts allows evaluating the ability of therapeutics on rescuing not only primary, but also secondary and tertiary neuro-pathological signatures.




Schematic of the AD model using iPSCs derived neurons, astrocytes and microglia treated with TAU PFFs. Day 10 post thaw co-cultures expressing 4R TAU P301L were treated with TAU PFFs. On day 18, Ncyte microglia were added. Co-culture of neurons, astrocytes and microglia was kept for 21 days for endpoint analysis of pTAU, MC-1 positive TAU in neurons (by ICC-IF).

1. AD related triggers in microglia monocultures




A) Quantification of IL-1 β , IL-6, TNF- α and IL-10 in monoculture of microglia after stimulation with different concentrations of Tau PFFs or LPS (Assay Pos Ctrl) or Cytochalasin D (Assay Neg Ctrl). Microglia is cultured in cell specific medium (white bars) or in tri-culture medium (purple bars). Median $\pm$ SD. Significance to respective vehicle control. **** p. < 0.0001

2. Ncyte® Neural Mix – Mature neurons and astrocytes co-culture expressing SYN and PSD-95




3. Establishment of CNS tri-culture system




Under basal condition, Ncyte microglia do not secrete key pro-inflammatory cytokines either in mono- or tri-culture systems. Upon stimulation with LPS, Ncyte® Microglia exhibit a strong pro-inflammatory reaction both in monoculture as well as in a tri-culture system comprising of neurons, astrocytes, and microglia. LPS treatment significantly increases the production of the key inflammatory mediators IL-6, IL-10 and TNF- α in both models, as measured by Mesoscale Discovery.

4. Establishment of tauopathy model in tri-culture system



High content-based quantification of the normalized parameters (pTAU puncta counts) of transduced tricultures (in black) compared to transduced tricultures treated with TAU PFFs.

5. Cytokine release in tri-culture system



Quantification of IL-1 β , IL-6, TNF- α , and IL-10 in monocultures of microglia following stimulation with increasing concentrations of TAU PFF or LPS (assay positive control). Cytokine levels were measured in the culture supernatant after 24 h using multiplex immunoassays. Data are presented as median $\pm$ SD from n = 3 biological replicates. TAU PFF stimulation induced a dose-dependent increase in pro-inflammatory cytokines (IL-1 β , IL-6, TNF- α) along with elevated levels of the regulatory cytokine IL-10, indicative of an activated microglial phenotype. Significance was determined relative to respective vehicle controls (**** p < 0.001; **** p < 0.0001).

Treatment with Tau PFFs induces release of pro-inflammatory cytokines, typical of AD pathology

Conclusions

- Microglia cultures exposed to TAU PFFs species released higher levels of pro-inflammatory cytokines IL-6, IL-8, IL-1 β and IL-10, an AD-related biomarker.
- We observe increased levels of pTAU and MC-1⁺ TAU in a mature tri-culture, as evident from the co-localisation of pre and post synaptic markers.
- Microglia in a tri-culture system demonstrate a significantly different activation profile compared to microglia in mono-culture indicating the need to include more complex models into your drug discovery project
- Co-culturing neurons, astrocytes, and microglia enhances the detection of Tau-induced neurodegenerative phenotypes, offering a more physiologically relevant platform for assessing therapeutic efficacy.

Evaluate the efficacy of your therapeutic candidates in mitigating neuroinflammation and protein aggregation using Ncardia’s human in vitro tri-culture model.

Contact us at
support@ncardia.
com or scan the
QR code to
download the
poster

