



Abstract ID: 16
Poster: #105

Predicting which healthy tissues an antibody-drug conjugate will damage in patients remains one of the most consequential unsolved problems in ADC development. Here we show that a fully human, iPSC-derived seven-tissue panel can reproduce the known clinical toxicity signature of Cetuximab-MMAE — including tissue selectivity, dose-dependency, and the absence of toxicity in tissues where it is not clinically observed — without reliance on animal models or immortalized cell lines. Tissue selection was informed by *in silico* EGFR mRNA expression profiling across human tissues (GTEx database), yielding a panel of seven iPSC-derived cell types spanning skin (melanocytes), liver (hepatocytes), lung (alveolar epithelial cells), heart (cardiomyocytes), CNS (motor neurons), vasculature (endothelial cells), and intestine (epithelial cells).

- Uses exclusively human iPSC-derived cell types
- Covers seven physiologically relevant tissues
- Enables real-time functional toxicity assessment
- End-point is non-invasive, allowing further optional analysis (i.e. omics)
- Provides mechanistic insight into tissue-specific susceptibility
- Builds upon prior IND-supporting experience in cell therapy programs

Desired activity		No activity		Unwanted activity	
Cancer cells Antibody Drug Conjugate  • Tumor target antigen recognition • Antibody internalization • Payload mediated cytotoxicity	iPSC-derived cardiomyocytes Antibody Drug Conjugate  • No target antigen recognition • No off-target antigen recognition • No antibody internalization • No payload mediated cytotoxicity	iPSC-derived Cells Antibody Drug Conjugate  • Target antigen recognition • Off-target antigen recognition • Antibody internalization • Payload mediated cytotoxicity			
Working Product	Safe Product		Unsafe Product		

```

graph LR
    A[In silico analysis] --> B[Validate target expression]
    B --> C[Impedance based cytotoxicity testing]
    C --> D[Optional molecular or biomarker readouts]
    D --> E[IND enabling report]
  
```

In silico analysis

Validate target expression

Impedance based cytotoxicity testing

Optional molecular or biomarker readouts

IND enabling report

Cell growth

0 days
1 day
2 days
3 days

Cell differentiation

0 days
1 day
2 days
3 days

Cell types: Cardiomyocyte, Hepatocytes, Spinal Motor Neurons, Melanocytes, Endothelial, Alveolar epithelial, Intestinal epithelial

xCELLigence MP

iPSC derived Lung alveolar epithelial cells

iPSC derived Vascular Endothelial cells

ADC dose-dependent cell killing

No ADC induced cell killing

Legend:

- Pos control ADC 100 µg/ml
- Pos control ADC 31.6 µg/ml
- Pos control ADC 10 µg/ml
- Pos control ADC 3.2 µg/ml
- Pos control ADC 1 µg/ml
- Pos control ADC 0.32 µg/ml
- Control IgG 100 µg/ml
- Control IgG 31.6 µg/ml
- Control IgG 10 µg/ml
- Control IgG 3.2 µg/ml
- Control IgG 1 µg/ml
- Control IgG 0.32 µg/ml
- DPBS

Tissue	EGFR Expression	In Vitro Toxicity (100 µg/mL)	In Vitro Toxicity (1 µg/mL)	Clinical AE Signal
Skin (Melanocytes)	High	High	Low residual	High
Liver (Hepatocytes)	Moderate	Moderate		Low–Moderate
Lung (Alveolar epithelial)	Moderate	Moderate	None	Moderate
Heart (Cardiomyocytes)	Low	Minimal	None	Minimal
CNS (Motor neurons)	Low	Minimal	None	Minimal
Vascular (Endothelial)	Low	None	None	Minimal
Intestine (Epithelial)	Low	None	None	Minimal